

## MRSD Team C - Fall Validation Test:

- **Location:** NSH B512
- **Objective:** Demonstrate that the system is capable of autonomously localizing, planning, and executing an acetabular reaming operation as it would be utilized in an operating room.
- **Equipment:** Atrascys camera and Vesa mount, Sawbone pelvis and vise, Kinova Gen-3 robot arm, Vention table, reaming end-effector
- **Personnel:** One team member interacting with the robot arm and work environment, one team member monitoring the Surgeon I/O and presenting.
- **Requirements:** M.P.1.1, M.P.1.2.1, M.P.1.2.2, M.P.2, M.P.3.1, M.P.3.2, M.P.4.1, M.P.4.2, M.P.4.3, M.P.5, M.P.6, M.P.7, M.P.8
- **Procedure:**
  - a. Begin by setting up the work environment by clamping the Sawbone pelvis in a new position in a vise, fixing a fiducial marker screw mount on the pelvis, and placing the fiducial marker onto the end-effector of the robot arm.
  - b. Utilizing a probe, the pelvis will be localized using a point cloud to fit the pelvis to a known pelvis mesh, from which the end-point of the reaming operation will be determined
  - c. Utilizing free motion mode, the robot arm will be placed near the center of the acetabulum.
  - d. The reaming operation would then be started, allowing the robot arm to localize itself with respect to the pelvis and begin generating a motion plan.
  - e. Once the reaming motor turns on and the arm begins to move, making contact with the pelvis, the e-stop is hit to demonstrate the safety of the system.
  - f. The robot arm will then be reset with free motion mode and the reaming operation would then be allowed to progress freely.
  - g. As the robot arm begins to ream the acetabulum, the pelvis would be shifted by hand using the vise, to demonstrate the robot arm's capability of adapting to pelvic motion.
  - h. When the robot arm has completed the reaming operation, it will remove itself from the pelvis, and the resulting acetabulum can be analyzed.
- **Validation:**
  - Pelvis pose must be localized with a latency  $\leq 50$ ms.
  - End-effector pose must be localized with a latency  $\leq 50$ ms.
  - Error in pelvis pose must be calculated with a latency  $\leq 50$ ms, a position error of  $\leq 1$  mm and orientation error  $\leq 1.5$  degrees.
  - Personnel can move robot arm freely during free motion mode
  - Once the e-stop is pressed, the system stops moving and the reamer stops spinning within 500 ms
  - Trajectory is generated within 150 ms
  - New trajectory is planned when pelvis is moved more than 1 mm or more than 1.5 degrees, and the execution is using the new trajectory
  - Using the Quantitative Trajectory Evaluator, evaluate the desired trajectory (from the trajectory planner) with the actual trajectory (recorded by the camera and joint encoders). The error must be below 1 mm and 1.5 degrees
  - Surgeon I/O shows progress and performance metrics during the reaming