

April 27, 2022



HIPSTER

Spring Validation Demonstration - Encore

Autonomous Reaming for Total Hip Replacement
(ARTHUR)



The Team



Kaushik Balasundar

Perception and
Sensing Lead



Parker Hill

Mechanical
Systems
Engineering Lead



Anthony Kyu

Controls and
Actuation Lead



Gunjan Sethi

Software
Engineering Lead



Sundaram Seivur

Trajectory Planning
Lead

Contents

(what's changed)



- Use Case
- Project Overview
- *Workspace*
- Recorded Tests
 - Hand-Eye Calibration
- Live Tests
 - Free Motion Mode
 - *Sparse Pointcloud Collection*
 - Controls
 - *Force feedback for Turning on Reamer*
 - Dynamic Compensation
 - Reaming the Pelvis
- Discussions and Questions


Use Case

Of the 100 manual surgeries, **30-45% of them observed the implant within the Lewinnek safe zone** and of the 100 robotic-assisted surgeries, **77% were within the safe zone.**


01

Study on the future projections on the number of total joint replacements in the US, show that up until 2040, we can expect an **increase in the requirement of a THR for both sexes by approximately 280%.**

02



Manual Reaming Demonstration




Overview



A **fully autonomous** robotic arm aimed at performing acetabular reaming with **high accuracy**, eliminating the need of surgeons to use intuition to correctly position/angle the reamer.

Workspace






Improvements for Encore



Force
Feedback

Improved
Motor Mount



More realistic
surgical setting

Tests

- Hand-Eye Calibration (Recorded)
- Free Motion Mode
- Pointcloud Collection
- Landmark Selection + Registration
- Planning
- Controls
- Dynamic Compensation
- Reaming the Pelvis



Hand-Eye Calibration



- **Objective:** Find transformation between /base_link (world) frame of robot and /camera frame
- Eye-on-base problem
- Marker used as calibration target
- Calibration done using OpenCV library's Tsai-Lenz algorithm implementation

Activities
Apr 19 14:37

Form

Info

Name:

Type:

Tracking Base Frame:

Tracking Marker Frame:

Robot Base Frame:

Robot Effector Frame:

Calibration algorithm:

Actions

Take Sample

Remove Sample

Compute

Save

Samples

```

1)
hand->world
translation: [-0.87, +0.05, -0.41]
rotation: [-0.57, +0.44, -0.47, +0.51]
camera->marker
translation: [+0.39, -0.03, +1.05]
rotation: [-0.13, +0.04, +0.98, -0.09]

```

Result

visualization.rviz - RViz

File Panels Help

Interact Move Camera Select Focus Camera Measure 2D Pose Estimate 2D Nav Goal Publish Point

Time

ROS Time: 1650393425.54 ROS Elapsed: 82.46 Wall Time: 1650393425.57 Wall Elapsed: 82.37

Reset Left-Click: Rotate. Middle-Click: Move X/Y. Right-Click/Mouse Wheel: Zoom. Shift: More options.

Experimental 31 fps

Hand Eye Calibration / SVD April 20 2022

Free-Motion Mode

Free motion mode will allow the surgeon to **move the end-effector to the patient's acetabulum** before executing the trajectory. The robot arm stays in place unless the surgeon moves it by hand.





Of the 100 manual surgeries, 30-45% of them observed the implant within the Lauenroth safe zone

Of the 100 manual surgeries, 77% of them observed the implant within the Lauenroth safe zone

The patient moves due to reaming forces.

The robot automatically aligns with the moved pelvis (dynamic compensation)

Anatomical Terms


Acetabulum: A cup-shaped socket that forms the ball-and-socket hip joint with the femur. The location we are reaming for a cup implant. The acetabulum is formed by the ilium, ischium and pubis extending from the acetabulum. The acetabulum is fixed, and the femoral head is positioned and anteversion of $40 \pm 10^\circ$ and $15 \pm 10^\circ$ to reduce the chance of dislocation when fitting a cup. The acetabulum must meet our system's performance requirements. The acetabulum is a spherical cutting shell which can fit into an acetabulum of a desired shape for an acetabular cup implant.


HMPSTER

Pointcloud Collection






*registration probe
(for pointcloud collection)*





Landmark Selection + Registration







Latency & Error Detection


○ Latency Test Procedure

- i. Clamp test model to Vention table.
- ii. Place a fiducial marker on the robot's end-effector.
- iii. Record the end-effector marker's pose from the camera.
- iv. Record time to get end-effector pose.

○ Error Detection Test Procedure

- i. Place a marker in the initial slot on the test model. Record its pose.
- ii. Move the marker to the planar slot. Record its pose.
- iii. Move the marker to the slanted slot. Record its pose.
- iv. Record computed translation and orientation error for the new marker positions.

```
min: 0.010s max: 0.026s std dev: 0.00097s window: 379
average rate: 54.033
min: 0.010s max: 0.026s std dev: 0.00092s window: 433
average rate: 54.035
min: 0.010s max: 0.026s std dev: 0.00088s window: 487
average rate: 54.032
min: 0.010s max: 0.026s std dev: 0.00085s window: 541
average rate: 54.034
min: 0.010s max: 0.027s std dev: 0.00095s window: 595
average rate: 54.036
min: 0.010s max: 0.027s std dev: 0.00092s window: 649
```





Latency & Error Detection

Performance Requirement:

Position Error $\leq 3\text{mm}$

Orientation Error $\leq 3\text{ degrees}$



Latency $< 500\text{ ms}$

Current System Performance:


Position Error $\leq 2\text{mm}$

Orientation Error $\leq 3\text{ degrees}$

Latency $\sim 20\text{ ms}$



Planning, Controls and Reaming (Encore)





Reaming Operation

Dynamic Compensation




ARTHuR constantly **checks for any movement of the patient** above a certain threshold and **re-plans the trajectory** of reaming if that threshold were to be crossed.






Trajectory Evaluation



Test 1






Test 2

	Req	Current
x	3mm	1.5mm
y	3mm	2.4mm
z	3mm	1.8mm



Trajectory Evaluation



Orientations threshold: **3 degrees (rpy)**

Hardware Design






End-Effector Design

Updated Design:

- Shortened original design to increase robustness and accuracy
- Enclosed motor housing for increased sanitization

Parts Necessary for Design:

- Force-Torque sensor
- Motor which can run at 400 rpm while outputting 0.5 Nm of torque
- Fiducial marker mount for camera-base registration
- Reamer head
- Reamer head adapter
- Motor mount
- Force-Torque adapter






PCB Design

Functionality:

- Needed to be able to turn on and control the motor utilizing ROS
- Further, needed the motor to maintain a consistent rpm during the reaming operation

Parts:

- Pololu 1457 -> Cytron MD10C
- Arduino Nano
- Header Blocks
- Fuses
- LEDs
- Resistors
- Power Supply



Changes for Fall Validation



Changes for Fall Validation

- Full **User Interface** which allows surgeons to easily interface with the system
- Full **PCB redesign** to integrate all components into one package
- Redesign to the **end-effector**:
 - Want to try different motor orientations (perpendicular orientation)
 - Potentially have stabilization rods
 - Machine out of aluminum or manufacture from plastics
- **Dynamic compensation** which occurs during the reaming operation (no replanning)
- Improved overall performance:
 - Reaming sawbones (closer bone approximate)
 - Faster planning
 - More accurate execution
- If time permits, implement **model predictive control** into our control system



Thank You!
See you next
semester :)

Autonomous Reaming for Total Hip Replacement (ARTHuR)

